


Establishing a robust platform for investigating the pro- and anti-inflammatory response of iPSC-derived microglia in drug discovery



*S. Humphreys, D. Wallbank, J. Tilman, H. Sharplin, H. Raine, S. Barreto
Axol Bioscience Ltd. B260 Meditrina, Babraham Research Campus, Babraham CB22 3AT; email: operations@axolbio.com

Abstract

Microglia are the main inflammatory cell of the brain and have been implicated in the development and progression of several neuroinflammatory and neurodegenerative conditions such as Alzheimer's disease (AD), Frontotemporal Dementia (FTD) and Amyotrophic Lateral Sclerosis (ALS). iPSC-derived microglia provide a scalable, reproducible and physiological relevant model with which to study this role. Axol Bioscience has developed a robust protocol to generate microglia, with over 30 lines successfully differentiated, including those derived from healthy and patient donors and gene-edited lines.

Microglia are highly plastic cells influenced by environmental cues to drive both pro- and anti-inflammatory responses. As such, any dysregulation in their activation status can drive an exaggerated and sustained state of inflammation. axoCells™ iPSC-derived microglia provide a reliable platform to study phenotypes of healthy and diseased microglia and assess the ability of novel compounds to decrease inflammation. In this study, Axol used a range of assays to investigate the activation state of microglia in response to different stimuli, as well as assessing inter batch and inter assay variability of our control (ax0664) microglia.

HTRF assays for IL8, IL6 and TNFα release demonstrated consistent effects across batches of microglia, when stimulated with a concentration range of lipopolysaccharide (LPS). In addition, HTRF and ELISAs were used to assess activation states of microglia induced by IL-4. Further characterization by flow cytometry, investigated the effect of different stimuli and incubation times on CD80, CD83 and CD206 receptor expression. Additional profiling of ax0664 investigated the release of a broad range of cytokines, using O-link technology. Furthermore, Tempo-Seq analysis was used to investigate transcriptomic profiles driven by different stimuli. Overall, this supports the ability of Axol's human iPSC-derived microglia to polarise towards different activation states. axoCells iPSC-derived microglia can be used as a relevant platform for human drug discovery where they can help predict the action of new drugs and identify mechanisms involved in neuroinflammation.

Methods

Human iPSC lines were differentiated to monocyte-like cells by mesoderm induction via embryoid-body formation ⁽¹⁾. After maturation towards microglia-like cells, cryopreservation of the cells was performed.




Figure 1. Flow diagram of differentiation from iPSC through a haemangioblast generation to cryopreserved microglial-like cells.

Upon thaw, quality control via immunocytochemistry staining for three key markers was performed to determine the identity of the cells.

Once the cells had passed QC, the microglia were treated with pro- / anti-inflammatory stimuli before measuring the IL8, IL6, TNFα and MMP12 concentration in the cell culture medium via HTRF or proteomics. In addition, TNFα, IL6 and IL8 release was measured following incubation of the cells with TLR4 antagonists + LPS.

Transcriptomics analysis, via TempoSeq, was performed on stimulated microglia. Finally, Flow cytometry following stimulation of the microglia was also performed to investigate the expression of three cell surface markers, CD80, CD83 and CD206.

Results

Quality Control of iPSC Microglia





Figure 2. Immunocytochemistry for three markers, IBA1, TMEM119 and P2RY12, on ax0664 microglia-like cells after seven days of maturation post thaw.

Microglia Proteomics



Microglia Transcriptomics

