

Culture medium
for human sebocytes
PhenoCULT®-SEB

Product sheet

DESCRIPTION

Product Ref: PhenoCULT-SEB

Culture medium for human sebocytes. Use according to the PCi-SEB_User's Guide for thawing and culture of sebocytes derived from iPSC (PCi-SEB).

PhenoCULT®-SEB is available in 100 mL or 500 mL format and contains:

Product	Label	Quantity
Basal sebocyte culture medium Supplements A (amplification) and M (maturation) 100 μ L or 5x100 μ L	PhenoCULT®-SEB_Basal	100 mL or 500 mL

PhenoCULT®-SEB has been optimized for iPSC-derived sebocytes and human primary sebocytes.

Associated Products:

- PCi-SEB: sebocytes from human iPS cells (1M cells/vial), available in 3 phototypes (Asian, Caucasian, African)
- PCi-SEB_KIT: sebocytes from human iPS cells (1M cells/vial), available in 3 phototypes (Asian, Caucasian, African) and PhenoCULT®-SEB

STORAGE AND EXPIRATION

- Storage conditions: < - 20° C. Store upon receipt.
- Expiration: refer to expiration date for each lot. Use within 6 months.

PRODUCT USE

PhenoCULT®-SEB is intended for IN VITRO RESEARCH USE ONLY and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

SAFETY PRECAUTIONS

These recommendations are based on prudent application of professional judgment. Wear the appropriate personal protection equipment (PPE) and handle with due caution. Do not

ingest. Safety Data Sheets are available on demand.

Phenocell cannot be held liable for any damage or losses resulting from the handling or from contact with the product as described herein.

QUALITY CONTROL

Refer to lot-specific Certificate of Analysis. PhenoCULT®-SEB is tested for bacterial, fungal and mycoplasma contamination, and for pH.

FOR RESEARCH USE ONLY

Not intended for human or animal diagnostic, therapeutic or clinical applications.